



**Al-Hami**

Medical Equipment LLC

Products Catalogue

**Electrophoresis**



**Al-Hami**  
Medical Equipment L.L.C



The new language of life

**ELECTROPHORES**

Capillarys 3  
TERA



<b>Method</b>	<b>Capillary electrophoresis</b>		
<b>Test Menu</b>	Protein, Hemoglobins, Immunotyping, CDT & Hb A1c		
<b>No. of Capillaries</b>	12	<b>Sample Capacity</b>	120 Sample
<b>Detection</b>	UV Absorbance	<b>Throughput</b>	Up to 100T/H



**Al-Hami**  
Medical Equipment L.L.C



The new language of life

**ELECTROPHORES**

Capillarys 3  
OCTA



<b>Method</b>	<b>Capillary electrophoresis</b>		
<b>Test Menu</b>	Protein, Hemoglobins, Immunotyping, CDT & Hb A1c		
<b>No. of Capillaries</b>	8	<b>Sample Capacity</b>	120 Sample
<b>Detection</b>	UV Absorbance	<b>Throughput</b>	Up to 80 T/H



**Al-Hami**  
Medical Equipment L.L.C

**sebia**  
The new language of life

**ELECTROPHORES**

Capillarys 2  
Flex Piercing



<b>Method</b>	<b>Capillary electrophoresis</b>		
<b>Test Menu</b>	Protein, Hemoglobins, Immunotyping, CDT & Hb A1c		
<b>No. of Capillaries</b>	8	<b>Sample Capacity</b>	Up to 80
<b>Detection</b>	UV Absorbance	<b>Throughput</b>	80-100 T/H



**Al-Hami**  
Medical Equipment L.L.C



The new language of life

**ELECTROPHORESIS**

miniCap  
Flex Piercing




<b>Method</b>	<b>Capillary electrophoresis</b>		
<b>Test Menu</b>	Protein, Hemoglobins, Immunotyping, CDT & Hb A1c		
<b>No. of Capillaries</b>	2	<b>Sample Capacity</b>	36
<b>Detection</b>	UV Absorbance	<b>Throughput</b>	10-15 T/H



**Al-Hami**  
Medical Equipment L.L.C

**sebia**  
The new language of life

**ELECTROPHORESIS**



Capillarys 2

<b>Method</b>	<b>Capillary electrophoresis</b>		
<b>Test Menu</b>	Protein, Hemoglobins, Immunotyping, CDT		
<b>No. of Capillaries</b>	8	<b>Sample Capacity</b>	Up to 80
<b>Detection</b>	UV Absorbance	<b>Throughput</b>	80-100 T/H



**Al-Hami**  
Medical Equipment L.L.C



The new language of life

**ELECTROPHORESIS**



miniCap

<b>Method</b>	<b>Capillary electrophoresis</b>		
<b>Test Menu</b>	Protein, Hemoglobins, Immunotyping, CDT		
<b>No. of Capillaries</b>	2	<b>Sample Capacity</b>	36
<b>Detection</b>	UV Absorbance	<b>Throughput</b>	10-15 T/H



Al-Hami

Medical Equipment LLC

[www.alhamico.com](http://www.alhamico.com)



[info@alhamico.com](mailto:info@alhamico.com)



+971-50 - 8507636