

## Waters e2695

Pump Type	Quaternary Gradient Pump
Flow Rate Range	0.001 – 10.000 mL/min
Flow Accuracy	±1% or ±0.01 mL/min
Pressure Range	Up to 5000 psi (≈345 bar)
Gradient Formation	Low-pressure mixing
Degasser	Built-in vacuum degasser
Autosampler Type	Integrated autosampler
Injection Volume	0.1 – 100 µL
Sample Capacity	Up to 120 vials
Column Compartment	External (optional heater)
Detector Compatibility	UV, PDA, Fluorescence, RI
Common Detector	Waters 2489 UV/Vis Detector
Software	Empower / MassLynx



Al-Hami  
Medical Equipment LLC

